
Monitor - Benutzerhandbuch

MDRAP-Series

Model: MDRAP21Exxxxxxx
(RAP2123AMII)

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
RAP Series

Model: RAP2123AMII

Thank you for purchasing our color LCD Display.

- ◆Carefully read this User's Manual and use the product properly. Before using it, also read "Safety Precautions."
- ◆Keep the User's Manual as close to you as possible and in safe custody.
- ◆If you have lost the manual, contact your dealer. We will reissue a manual.
- ◆The names of companies and products are registered brand names or brand names.

Contents	
1.Safety Precautions.....	1
2.Equipment Symbols.....	4
3.Auxilianies.....	6
4.Application & Function Description.....	7
5.Features.....	8
6.Connection Method.....	9
7.How to Operate.....	11
8.Applicable Signals.....	19
9.Failure Diagnosis.....	21
10.Cleaning.....	22
11.Specifications.....	23
12.Technical Support.....	24
13.Contact Information.....	25

Revision History

The revision letter changes with related comments each time the document is updated.

Revision	Comment
A	Initial release of this document

Safety Precautions

I CAUTION

Symbol Description

Caution!-This symbol alerts you to important operating considerations or a potential operating condition that could damage equipment. Refer to user's manual or operation's manual for precautionary instructions.

II CLASSIFICATION

MEDICAL ELECTRICAL EQUIPMENT
WITH RESPECT TO ELECTRICAL SHOCK, FIRE AND
MECHANICAL HAZARDS ONLY IN ACCORDANCE
WITH ANSI/AAMI ES60601-1 (2005,3rd ed.) AND
CAN/CSA C22.2 No. 60601-1 (2008,3rd ed.)
Classification: Class I
No Applied parts
No AP/APG

Class I :

No applied parts

Protection against harmful ingress of water is IPX0

Not suitable for use in the presence of flammable anesthetic's or oxygen.

Mode of operation: Continuous.

III External Equipment

External equipment intended for connection to signal input / output or other connectors, shall comply with UL/EN 60601-1 for medical electrical equipment. In addition, all such combination system shall comply with the standard IEC 60601-1-1, Safety requirements for medical electrical systems. Equipment not complying with UL/EN/IEC 60601-1 shall be kept outside the patient environment, as defined on the systems standard.

IV Intended Use

The equipment is intended to be used as a component of a medical patient monitoring system.

V FCC

FCC Information

FCC(U. S. Federal Communications Commission)

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses, and can radiate radio frequency energy, and if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause unacceptable interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult your dealer or an experienced radio/TV technician for help.

FCC Warning:

To assure continued FCC compliance, the user must use a grounded power supply cord and the provided shielded video interface cable with bonded ferrite cores. Also, any unauthorized changes or modifications to this monitor would void the user's authority to operate this device.

VI Precautions

- Make sure to carefully read the Operating Manual prior to using the monitor in order to properly use the monitor.
- Note that, excluding those cases where a responsibility for legal compensation is recognized, the manufacturer shall bear absolutely no responsibility for damage to this product by a customer or a third party that results from the ignoring of contents entered in this Operating Manual and mistaken use.
- **Follow the instructions below for safety use of the LCD Monitor.**
 - To avoid electric shock, do not attempt to remove any cover or touch the inside of the monitor. Only a qualified service technician should open the monitor case.
 - Do not insert metal objects or spill liquid into the LCD monitor through cabinet slots. They may cause accident fire, electric shock or failure. If a foreign object inserted or water penetrated, unplug the AC cable and have the monitor serviced.
 - Do not cover or block the vent holes in the case.
 - Disconnect the power plug from the AC outlet if you will not use it for an indefinite period of time.
 - Do not touch the screen directly with your fingers. You may damage the screen, and oil from your skin is difficult to move.
 - Do not apply pressure to the screen. The LCD is very delicate.

- If your LCD monitor uses an AC/DC Adapter: Only use the Adapter, which accompanied with this device. Use of another AC/DC Adapter may cause a malfunction or danger.
- If your LCD monitor does not operate normally-in particular, if there are any unusual sounds or smells coming from it-**unplug the AC cable** immediately and contact the manufacturer, authorized or service center.

VII Environment

- Place the monitor on a flat and leveled surface.
- Place the monitor in a well-ventilated place.
- Keep the monitor away from:
overly hot, cold or humid places, places directly under sunlight, dusty surroundings, equipment that generate strong magnetic fields.

Equipment Symbols

Electrical and Electronic Equipment Symbols

In addition to the equipment symbols described in your user's manual, the following symbols may be appear on the monitor.

 Alternating current.

 Direct current.

 European Union Declaration of Conformity.

 FCC. USA only. Complies with applicable US government (Federal Communications Commission) radio-frequency interference regulations.

 Fragile. Handle with care.

 Indicates front.

 Keep dry. Protect from rain.

 Mercury. This product consists of devices that may contain mercury, which must be recycled or disposed of in accordance with local, state, or country laws.


 ON. Power connection to the mains.

 Recycled materials or may be recycled.

 Stacking limit by number.

 Standby or power indicator.

 This way up.

 USB connector port.

This symbol indicates that the waste of electrical and electronic equipment must not be disposed as unsorted municipal waste and must be collected separately. Please contact an authorized representative of the manufacturer for information concerning the decommissioning of your equipment

ATTENTION: Consult accompanying documents.

DANGER - Shock Hazard. Dangerous voltage. To reduce the risk of electric shock, do not remove cover. Refer servicing to qualified service personnel.

Equipotentiality. Connect device to a potential equalization conductor.

Auxiliaries

Check List

Before operating this monitor, please make sure that all items listed below are present in your package:

- The LCD monitor
- This user's manual
- AC Power cord(USA/GERM)
- AC /DC adapter
- D- sub 15 -pin signal cable
- DVI/DP cable
- Audio Cable

Optional Items

These are optional items and can be ordered by contacting dealer directly.

- VESA® standard wall/arm mount kits

- *1 USB / RS232 Cable

If any of the above items are found missing or if you wish to order the optional items, please contact your dealer.

*1 For the Touch unit with the USB / RS232 input connector

NOTES

* We recommend you to keep the packing box for transportation and the like.

* In case of transporting and packing the unit in the packing box, carefully place it, keeping its panel from touching any objects.

Application and Function Description

●Some Basic Facts

The medical safety approved RAP model is a versatile monitor designed for ergonomic and easy handling. The liquid-crystal display (LCD) monitor supports the most common resolutions from 640x480 (VGA) up to 1920x1080(HD 1080) and presents sharp, low-radiation images.

With the OSD wheel knob on the bottom right side of the display leaves ample space on your desk for other peripheral equipment.

An VESA 100 adapter according to VESA standard is available on the rear cover for wall or swivel-arm mounting.

●Function Description

The medical flat-screen display is built into a slim and ergonomic housing. The monitor displays 16.7 Million colors with a resolution of 1920x1080 pixels.

The signal is equivalent to the analog and digital standard signals of your PC.

After connection of the signal the flat-screen display automatically adapts, as far as possible, to the VGA signal and presents a stable and centered image. Additional display parameters, such as brightness and contrast, can be adjusted via the on-screen menu.

The equipment supplied includes the signal cables, the medical-grade power supply with power cord and this user manual.

Features

●21.5-inch LCD display with 2.1 million pixels

This color LCD display has a multi-scanning function corresponding to the resolution from VGA 640x400 to HD 1080 1920x1080. This is also compliant with VESA standard display mode.

●Power management

This unit has loaded the power management system. The power management mode functions when either horizontal or vertical signals or both disappears, and it reduces power consumption to less than 5W.

●VESA® standard wall / arm mountings

The unit is compliant with VESA's 100mm-pitch hanging tools. The tilt stand is detachable; the unit can be set for wall-hanging or arm according to users' environment.

●Security device slot for Kensington® security system

Anti-theft security device slot has been loaded to correspond to MicroSave® security system of Kensington®.

Note 1) PACS: Picture Archive & Communication System

Connection Method

No tools are required to set up the LCD monitor. Simply follow the instructions outlined in the next few pages.

Connectors for the signal cables, power are located on the back of the panel behind the cover door. Please refer to the diagram on this page for the connector configuration.

Connect Signal Cable(VGA)

Attach the VGA signal cable end, which has the ferrite core furthest away from the plug, to the monitor and attach the other end to the graphics card adapter on your computer. Be cautious in inserting the cable properly into both connectors. If the cable does not seem to fit it may be facing the wrong direction. Turn the cable over and try to match the shape of the connector with that of the graphics adapter.

Connect Power Adapter and Cable

Connect the round shape plug end of the AC/DC adapter to the DC Power input connector of the LCD monitor. Connect the female end of the power cable to the AC power input receptacle on the AC/DC adapter. Then, plug the male end of the power cable into an AC outlet.

Connect DVI/DP Cable

If you have a DVI/DP digital graphics card adapter and a DVI/DP cable, connect it to the DVI/DP connector of the monitor. The optional DVI/DP graphics adapter and the DVI/DP cable can be ordered by contacting dealer directly.

Connect Audio Cable

Connect Audio in connector of the LCD monitor and Audio out device (PC,DVD,CD...) with the attached stereo mini cable.

CAUTION: Make sure you use the auxiliary AC adapter.

***1 Connect USB / RS232 Cable**

Connect USB / RS232 in connector of the LCD monitor and PC out device. With the attached USB / RS232 cable.

CAUTION: When you disconnect the cord /cables, be sure to hold the connector and not the cable itself.



Cable Management

The LCD monitor has built-in cable guides and cover to help you organize and route the cables neatly on your desktop space. The cable cover is located on the back of the monitor stand and it is opened by gently pulling it off the stand. Under the cover, the cables can be held in position by numerous clips.

There are several ways to route the cables and you may use any way you find suitable for your situation. However, the routing method illustrated by the drawing on the right is recommended to give the best result in most situations both aesthetically and functionally. Once the cables are in position, gently snap the cover back on to the stand.

Anti-Theft Security Lock Slot

The LCD monitor is equipped with a security lock slot compatible to Kensington® security lock type. The security cable lock maybe available thru your dealer or it can be purchased at most computer peripheral stores near you.


Wall/arm VESA mounting the Monitor If you intend to mount the monitor on the wall/arm, we strongly recommend that you use only UL approved wall/arm mount kits with attached M4*10mm screws and can load more than monitor weight, that you ensure it is securely and safely installed. (VESA mounting for 100*100mm).

If you use a non-UL approved wall/arm mount kits, there is a safety risk that the Monitor may fall from the wall/arm due to an improper attachment through the use wrong or length screws.

How to Operate

Names and Functions of Each Part

Control Dial

The Control Dial is a multi-functional device located behind the LED Indicator on the right side of the front bezel. It has three movements-rotate upward, rotate downward and press inward as a button.

1.Power On/Off

Press the Control Dial to power the unit on from the off stage (the LED is off). To turn the power off, press the Control Dial and hold for at least 1 second until the LED turns off.

2.OSD Control

While the monitor is on (green LED and image on the screen), Pressing on the Control Dial activates the OSD. While the OSD menu is active, use the three way movements of the Control Dial to adjust the monitor.

Rotate Downward: Move Up/Right, Increase, Larger, More

Rotate Upward: Move Down/Left, Decrease, Smaller, Less

Button Press: Execute, Do, Save

LED Indicator

This LED indicator turns green when the power is switched ON and the power cord is properly attached. It turns amber when the monitor goes into a power saving mode (Active Off). Please refer to the Power Management section of manual for more information.

Auto Setup*

Automatically adjust the size, position, contrast, and the like of the screen.

Is this image displayed properly?

Yes: Return to Main Menu

No: Return to Clock/Phase (Change manually)

Brightness

Adjust brightness of the screen. The larger the value is, the screen is , and vice versa.

Adjustable at the range from 0 to 100 by a couple of steps.

Contrast

Adjust contrast of the screen. The larger the value is, the screen is , and vice versa.


Adjustable at the range from 0 to 100 by a couple of steps.

Display*


Exit

Return to Main Menu.


H. Position*

Adjust the horizontal position of the display. When rotating control dial clockwise, the display moves to the right, and vice versa. For the digital input DVI or Display Port, this function is not available.

V. Position*


Adjust the vertical position of the display. When rotating control dial clockwise, the display moves to the up, and vice versa. For the digital input DVI or Display Port, this function is not available.


Color Mode

User can make tone adjustments and settings.

- **Cool:** Bluish white for general use and CAD / CAM.
9300K is a fixed value by the manufacturer; User cannot change the setting value.
- **Neutral:** White close to natural light used for publishing trade.
6500K is a fixed value by the manufacturer; User cannot change the setting value.
- **Warm:** Reddish white mainly suited for photo modification. 5400K is a fixed value by the manufacturer; User cannot change the setting value.


- **User:** User can adjust and set tones.
 - **Red:** Adjust red and equivalent colors at the range from 0 to 100.
 - **Green:** Adjust green and equivalent colors at the range from 0 to 100.
 - **Blue:** Adjust blue and equivalent colors at the range from 0 to 100.

Clock / Phase*


Exit

Return to Main Menu.


Clock*

When operating "Auto adjustment" and if any partially phased screen, make fine adjustments on video signal position to horizontally synchronized signals by rotating the control knob. As the value gets larger, the screen moves to the right, and vice versa. Adjustable at the range from 0 to 100. For the digital input DVI or Display Port, this function is not available.


Phase*

Make "Phase " adjustments, if any flickers, blurs, or horizontal stripes on screen. Adjustable at the range from 0 to 63. For the digital input DVI or Display Port, this function is not available.


Management

OSD Display


Exit

Return to Main Menu.


OSD H. Position

Adjust the horizontal position of OSD. When rotating control dial clockwise, OSD moves to the right, and vice versa.

OSD V. Position

Adjust the vertical position of OSD. When rotating control dial clockwise, OSD moves to the up, and vice versa.

Language

Selecting this control to select the language you want.

- **English:** Display OSD in English.
- **Francais:** Display OSD in Francais.
- **Deutsch:** Display OSD in Deutsch.
- **Espanol:** Display OSD in Espanol.
- **Italiano:** Display OSD in Italiano.
- **こホン:** Display OSD in Japanese.


* Used in analog source input.

Source

Selecting this control allows you to manually select VGA or DVI or Display Port to automatically detect input sources all.

- **Exit:** Return to Main Menu.
- **VGA:** Display signal from VGA connector (15-pin Mini D-sub).
In VGA input, it supports Sync-on green and the interlace timing can be displayed.
- **DVI:** Display signal from DVI connector. (24-pin DVI)
- **Displayport:** Display signal from Displayport connector.

Recall


Initialize the data such as display position and automatic adjustment data to factory default. After initialization, perform "Auto Adjustment" again.

Mute

To control the audio sound of monitor.

ON: mute


Off: no mute

Volume

Adjust the volume of the audio sound.

Selecting this control to volume a scaling. At the range from 0 to 100.

Applicable Signals

1.Applicable Signals Timings

※The display may not work correctly with timings other than listed below.


◎ : Recommended timing

○ : Applied timing

Mode Name	Resolution		Frequency		Input Signal		
	Horizontal	Vertical	Horizontal (kHz)	Vertical (Hz)	Analog	Digital	DP
VGA	720	400	31.47	70	○	○	○
	640	480	31.47	60	○	○	○
VESA	640	480	37.86	72			
	640	480	37.50	75			
	800	480	35.16	56			
	800	600	37.88	60	○	○	○
	800	600	48.08	72			
	800	600	46.88	75			
	1024	768	48.36	60	○	○	○
	1024	768	56.48	70			
	1024	768	60.02	75			
	1280	1024	63.98	60	○	○	○
	1280	1024	79.98	75			
	1600	1200	75.00	60	○	○	○
	1920	1080	67.50	60	◎	◎	◎


2.Assignment of Connector Pins

15-pin Mini D-sub Connector (female)

VGA signal connector


B-Type Connector (female)

USB signal connector


Pin#	Signal	Pin#	Signal
1	S1(R:Red)	9	-
2	S2(G:Green)	10	HS-GND
3	S3(B:Blue)	11	VS-GND
4	GND	12	DDC,SDT
5	GND	13	HS
6	S1-GND	14	VS
7	S2-GND	15	DDC,SCK
8	S3-GND		

Pin#	Signal
1	VCC
2	D-
3	D+
4	GND

DIGITAL Signal Connector

Pin#	Signal
1	T.M.D.S. Data2-
2	T.M.D.S. Data2+
3	T.M.D.S. Data2/4 Shield
4	NC
5	NC
6	DDC Clock
7	DDC Data
8	NC
9	T.M.D.S. Data1-
10	T.M.D.S. Data1+
11	T.M.D.S. Data1/3 Shield
12	NC
13	NC
14	+5V Power
15	Ground
16	Hot Plug Detect
17	T.M.D.S. Data0-
18	T.M.D.S. Data0+
19	T.M.D.S. Data0/5 Shield
20	NC
21	NC
22	T.M.D.S. Clock Shield
23	T.M.D.S. Clock +
24	T.M.D.S. Clock -

RS232 signal connector

Pin#	Signal
1	NC
2	TXD
3	RXD
4	NC
5	GND
6	NC
7	NC
8	NC
9	NC

Failure Diagnosis

Symptom	Verify	Countermeasure
Power doesn't turn on! (POWER indicator doesn't light!).	a: Is the power cord connected? b: Is the power ON?	a: Connect the power cord(see pg.10) b: Turn the power on.
The screen is not displayed!*(POWER indicator lights up).	a: Is the system on? b: Is contrast or brightness adjustment at minimum? c: Is the signal cable properly connected? d: Is the signal input an applicable one within the stipulated frequency range?	a: Turn the system on. b: Adjust contrast and brightness(See pg.12). c: Properly connect the appropriate signal cable (see pg.9). d: Input the applicable signals within the frequency range.
Display size and position are not appropriate!	a: Is the preset adjusted to comply with input signals?	a: Make display adjustments (See pg.13~18).
The display vibrates!	a: Is the signal that is out of the set frequency range input?	a: Input applicable signals within the frequency range.
Control dial doesn't work!	a: Are you not rotating the dial so fast?	a: Rotate the dial slowly.

Note

*:Signals from computers may not come into when power management function is on.

Cleaning the Cabinet and the Color LCD Display

- When cleaning, remove the outlet from the LCD display for safety.
- Lightly wipe off dirt on the cabinet and LCD panel surface with a hard-wrung soft cloth soaked in a neutral cleaning solution. Follow its instruction when using a disposable cloth.
- Do not use thinner, benzine, alcohol or such on the cabinet that is made of plastic. These can damage the cabinet, alter its quality and cause the paint to peel off.
- Do not apply insecticides and other volatile items to the cabinet. Also do not leave rubber and vinyl products or such in contact with it for long hours. This can cause the quality to alter and the paint to peel off.
- Cleaners usable for the LCD panel are isopropyl alcohol(without abrasive), non-ammoniac glass cleaner, and watered-down neutral cleaning solution. Do not use organic solvent such as acetone and toluene.
- When the screen has dust on the LCD panel surface, wipe it off with soft moist cloth.
- Treat the LCD panel with care. Do not rub the LCD panel surface with a rough item or hit it on the surface. Also, do not strongly press the LCD panel surface. This can lead to unevenness in the screen and also to failure of the product.


Specifications

Product name :21.5"Color LCD Display

Items	Specifications		
LCD display device	55cm (21.5inch)Color TFT Normally Black		
Pixel pitch	Horizontal 0.248mm x Vertical 0.248mm		
Display area	Horizontal 476.6mm x Vertical 268.1mm		
Pixel	1920 x 1080 pixels		
Display gradation	16.77 million (8 bit each)x3		
Standard viewing angle	Horizontal : 178 deg. Vertical : 178 deg.		
Input signal	(1)VGA (15-pin Mini D-sub)connector.:0.7Vpp/75 Ohm, separate TTL level (2)DVI(24-pin) connector (3)Display port connector		
Temperature		On	Off
	Temperature	: 10~40℃	-20~60℃
	Humidity(non-condensation)	: 30~75%	10~90%
	Air pressure	: 700~1060hPa	500~1060hPa
Power Supply	DC 24V, 2.3A; (Provided Adapter: Delta Elec., Type: MDS-060AAS24 B, I/P: 100-250Vac, 1.5A-0.75A, 50-60Hz, O/P:24Vdc, 2.5A)		
Consumption current	Approx25W Max./Less than 2.5W when power management on (without Audio) Approx30W Max./Less than 2.5W when power management on (with Audio)		
External dimensions	Width 537.37mm x Depth 59mm x Height 328.87mm(landscape)		
Mass	Approx 5.7kg (net weight)		
International standards	UL/CUL (ANSI/AAMI ES 60601-1 and CAN/CSA C22.2 No.60601-1) CE, FCC-B		

External Dimensions

Landscape(longer side at the top)

●The specifications and appearance will be changed for improvement without notice.

Technical Support

Requesting Repair

- ① Read "Failure Diagnosis" carefully and check them yourself.
- ② In case of having abnormalities on your driver, stop operating, unplug the AC cable from the outlet, and Consult dealer. Do not repair by yourself. It is very dangerous.
- ③ Consult dealer for repair.
- ④ Use the product's packing and packing materials to send the product. Carefully place the unit in the packing, and keep the face from touching any packing materials.
The panel surface touching the materials during transportation will cause damage on the panel surface; we will not cover the damage.

Contact Technical Support

Consult your dealer for technical support.

Notes for User's Manual

- It is prohibited that copying any part or all of this manual without authorization.
- The content of this manual will be subject to change without notice.
- This was made carefully. Please contact us if any unclear points, mistakes, or parts left out.

Technical Support

Richardson Electronics GmbH

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D-78166 Donaueschingen

Telefon: + 4 9 - (0 7 7 1) – 8 3 0 0 – 4 7 / 4 9 .

Fax: + 4 9 - (0 7 7 1) – 8 3 0 0 – 8 0 .

www.canvys.de

Consult your dealer for technical support.

RAP LCD Display
Model : RAP2123AMII
ZCDDLE3010